



Product Datasheet

NME2 Antibody (orb229797)


Description	Mouse monoclonal antibody to NME2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Purified recombinant fragment of human NME2 (AA: FULL(1-152)) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	17.3kDa
Entrez	4831
Dilution Range	IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black line:
Control Antigen (100
n...



Immunohistochemical
analysis of formalin...



Western blot analysis
of human
recombina...