

Product Datasheet

phospho-PI3 Kinase p110 beta (Ser1070) antibody (orb2298)

Description phospho-PI3 Kinase p110 beta (Ser1070) antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IHC-P, WB

Immunogen KLH conjugated synthesised phosphopeptide derived from human PI3 Kinase p110 beta around the phosphorylation site of Ser1070 (YR(p-S))

Target phospho-PI3 Kinase p110 beta (Ser1070)

Preservatives 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance Liquid

Concentration 1mg/ml

Storage Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 110-123 kDa

Uniprot ID [P42338](#)

Dilution Range WB=1:500-2000, IHC-P=1:100-500, ICC=1:50

Expiration Date 12 months from date of receipt.



Western blot: Lane 1: Human MCF-7 cell; Lan...