



## Product Datasheet

### HSF1 Antibody (orb229806)


|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to HSF1                                               |
| <b>Species/Host</b>        | Mouse                                                                           |
| <b>Reactivity</b>          | Human                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Tested Applications</b> | ELISA, ICC, IHC                                                                 |
| <b>Immunogen</b>           | Purified recombinant fragment of human HSF1 (AA: 256-359) expressed in E. Coli. |
| <b>Preservatives</b>       | Purified antibody in PBS with 0.05% sodium azide.                               |
| <b>Concentration</b>       | 1mg/ml                                                                          |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.                     |
| <b>Note</b>                | For research use only                                                           |
| <b>Isotype</b>             | Mouse IgG2b                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                      |
| <b>MW</b>                  | 57.3kDa                                                                         |
| <b>Entrez</b>              | <a href="#">3297</a>                                                            |
| <b>Dilution Range</b>      | IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |



ELISA Black line:  
Control Antigen (100  
n...



Immunofluorescence  
analysis of A431  
cell...



Immunohistochemical  
analysis of formalin...