



Product Datasheet

MLANA Antibody (orb229807)


Description	Mouse monoclonal antibody to MLANA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human MLANA (AA: 48-118) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	13.2kDa
Entrez	2315
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black
line: Control
Antigen (100
n...



Western blot
analysis of
human
recombina...



Western blot
analysis of
HEK293 using
ML...