



## Product Datasheet

### PIK3CG Antibody (orb229811)


|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to PIK3CG                                             |
| <b>Species/Host</b>        | Mouse                                                                           |
| <b>Reactivity</b>          | Human                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Tested Applications</b> | ELISA                                                                           |
| <b>Immunogen</b>           | Purified recombinant fragment of human PIK3CG (AA: 1-100) expressed in E. Coli. |
| <b>Preservatives</b>       | Purified antibody in PBS with 0.05% sodium azide.                               |
| <b>Concentration</b>       | 1mg/ml                                                                          |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.                     |
| <b>Note</b>                | For research use only                                                           |
| <b>Isotype</b>             | Mouse IgG1                                                                      |
| <b>Clonality</b>           | Monoclonal                                                                      |
| <b>MW</b>                  | 126.5kDa                                                                        |
| <b>Entrez</b>              | <a href="#">5294</a>                                                            |
| <b>Dilution Range</b>      | ELISA: 1/10000                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |



ELISA Black  
line: Control  
Antigen (100  
n...



Western blot  
analysis of  
human  
recombina...



Western blot  
analysis of  
HEK293 using  
PI...