



Product Datasheet

RANGAP1 Antibody (orb229815)


Description	Mouse monoclonal antibody to RANGAP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human RANGAP1 (AA: 359-587) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	63.5kDa
Entrez	5905
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black line: Control
Antigen (100 n...



Western blot analysis of human recombin...



Western blot analysis of HEK293 using RA...