



Product Datasheet

Integrin alpha 4 antibody (orb229830)


Description	Integrin alpha 4 antibody
Species/Host	Rabbit
Reactivity	Bovine, Equine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Integrin alpha 4. Please contact us for the exact immunogen sequence. The peptide is available as orb374845.
Target	Integrin alpha 4
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	150 kDa
Uniprot ID	Q6PJE7, Q00651, D3ZMQ3
Dilution Range	WB: 2 µg/ml, IHC-P:1:100, IF/ICC: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of human breast ca...



Western blot analysis of rat uterus tiss...



Immunohistochemical staining of paraffin...