



Product Datasheet

CD90 antibody (orb229832)


Description	CD90 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CD90. Please contact us for the exact immunogen sequence. The peptide is available as orb374844.
Target	CD90
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P04216, P01831, P01830
NCBI	006288, 21, 31, 006279
Entrez	7070
Dilution Range	IF/ICC: 1:100, IHC-P:1:100
Expiration Date	12 months from date of receipt.



IHC-P image of human breast cancer tissue...



IF analysis of human breast cancer tissue...



IHC-P staining of human breast cancer tissue...