



Product Datasheet

CDK2 antibody (orb23096)


Description	CDK2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Cyclin-dependent kinase 2 protein (1-298AA)
Target	CDK2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24941
Dilution Range	WB:1:1000-5000, IHC-P:1:500-1000, IF:1:50-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562(lane 1), H...



Immunohistochemical staining of human to...



Immunohistochemical staining of human CO...