



Product Datasheet

SRRM2 antibody (orb2310)

Description	SRRM2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human SRRM2 (101-230/2752 aa)
Target	SRRM2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	299 kDa
Uniprot ID	Q9UQ35
Expiration Date	12 months from date of receipt.



Paraffin-embedded mouse liver tissue fix...



WB analysis of mouse, rat and human brea...