


Product Datasheet

CHM antibody (orb23175)

Description	CHM antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Rab proteins geranylgeranyltransferase component A 1 protein (1-653AA)
Target	CHM
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24386
Dilution Range	IHC-P:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...