


Product Datasheet

ABCG2 antibody (orb233596)

Description	Goat polyclonal antibody to ABCG2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ABCG2 / BCRP (Restricted sale due to demands of patent holder)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in lysates of cell line CACO-2. Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	EPSKQDKPLIEK
Clonality	Polyclonal
MW	72.3; 67.5
Entrez	9429
Dilution Range	ELISA: 1:1000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of CACO-2 lysate u...