


Product Datasheet

ALDOA antibody (orb233601)

Description	Goat polyclonal antibody to ALDOA
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ALDOA (aa86-96)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoformsWB: Approx 39kDa band observed in Human, Mouse and Rat Skeletal Muscle lysates. Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	QKADDGRPFQ
Clonality	Polyclonal
MW	39.4; 45.3
Entrez	226
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human, Mouse an...