



Product Datasheet

ASIC1 antibody (orb233606)


Description	Goat polyclonal antibody to ASIC1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ASIC1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all reported isoformsWB: Approx 70kDa band observed in Mouse fetal Brain lysates. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	NILPHHPARGT
Clonality	Polyclonal
MW	64.8; 59.9; 62.7
Entrez	41
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Co...



Western blot analysis of Mouse fetal Bra...



Immunohistochemical staining of Human Pa...