



Product Datasheet

ATP5C1 antibody (orb233609)

Description	Goat polyclonal antibody to ATP5C1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ATP5C1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoformsWB: Approx 33kDa band observed in Human Heart and in Rat Heart lysates. Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	TLTFNRTRQAVITKE
Clonality	Polyclonal
MW	32.9; 33.0
Entrez	509
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Heart lysat...



Immunohistochemical staining of Human He...