



Product Datasheet

ATP5F1 antibody (orb233610)

Description	Goat polyclonal antibody to ATP5F1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ATP5F1 (aa142-153)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 25kDa band observed in Human Heart lysates
Protein Sequence	SIQHIQNADTE
Clonality	Polyclonal
MW	28.9
Entrez	515
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...



Immunohistochemical staining of Human He...