


Product Datasheet

BMP1 antibody (orb233614)


Description	Goat polyclonal antibody to BMP1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	BMP1 (aa284-296)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_001190.1; NP_006120.1).ELISA: antibody detection limit dilution 1:128000.WB: Approx 110kDa band observed in Human Heart and Human Kidney lysates. Recommended concentration: 1-3 µg/ml.
Protein Sequence	DTIVPKYEVNGVK
Clonality	Polyclonal
MW	111
Entrez	649
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of rat Heart lysat...



Western blot analysis of NIH3T3 cell lys...



Western blot analysis of NIH3T3 cell lys...