


Product Datasheet

HS3ST1 antibody (orb233626)

Description	Goat polyclonal antibody to HS3ST1
Species/Host	Goat
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HS3ST1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 35kDa band observed in Human Brain lysates. Recommended concentration: 1-3 µg/ml.
Protein Sequence	NKLHEYFHEPNKK
Clonality	Polyclonal
MW	35.8
Entrez	9957
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...