


Product Datasheet

ITGB1 antibody (orb233627)

Description	Goat polyclonal antibody to ITGB1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ITGB1 / CD29
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 150kDa band observed in Human Gastrointestinal cancer and Gall Bladder cancer lysates. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	NKCPKKDSDSKIRP
Clonality	Polyclonal
MW	88.4; 87.4; 91.6; 91.0; 88.9
Entrez	3688
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Gastrointestinal...