


Product Datasheet

Nlrp12 antibody (orb233634)

Description	Goat polyclonal antibody to Nlrp12
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Nlrp12 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 110kDa band observed in Mouse Liver lysates. Recommended concentration: 1-3 µg/ml.
Protein Sequence	DDPPEPSGVQTQST
Clonality	Polyclonal
MW	119
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Liver lys...