


Product Datasheet

Sp1 antibody (orb233641)

Description	Goat polyclonal antibody to Sp1
Species/Host	Goat
Reactivity	Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Sp1(mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100kDa band observed in nuclear lysates of cell line NIH3T3
Protein Sequence	DGIHDETAGGEGEN
Clonality	Polyclonal
MW	80.4
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 nuclear ...