


Product Datasheet

TERF1 antibody (orb233643)

Description	Goat polyclonal antibody to TERF1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TERF1 / TRF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 65kDa band observed in Human Testis lysates. The observed aberrant MW may be explained by the acetylated form of the protein. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	NLDTRKSVSDKQS
Clonality	Polyclonal
MW	50.2; 48.2
Entrez	7013
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Testis ly...