



Product Datasheet

Vimentin antibody (orb233645)


Description	Goat polyclonal antibody to Vimentin
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	Vimentin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in lysates of cell line HeLa and Jurkat. Recommended concentration: 1-3 µg/ml.
Protein Sequence	QVINETSQHHDDLE
Clonality	Polyclonal
MW	53.7
Entrez	7431
Dilution Range	ELISA: 1:4000, WB: 0.1-2 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb233645
(2µg/ml)
staining of
HeLa (A)...



orb233645
(0.1µg/ml)
staining of
Mouse ...



Flow
cytometric
analysis of
paraformalde...