


Product Datasheet

FCRL2 antibody (orb233656)

Description	Goat polyclonal antibody to FCRL2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FCRL2 (aa92-104)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in lysates of Human Tonsil and of cell line Daudi. Recommended concentration: 1-3 µg/ml.
Protein Sequence	DKTSNIVKIKVQE
Clonality	Polyclonal
MW	55.5
Entrez	79368
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Tonsil (L...