



Product Datasheet

PLA2G2A antibody (orb233658)


Description	Goat polyclonal antibody to PLA2G2A
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PLA2G2A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Reported variants represent identical proteinWB: Approx 19kDa band observed in Human Placenta lysates. Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	SYKFSNSGSRIT
Clonality	Polyclonal
MW	16.1
Entrez	5320
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Placenta lysate...



Immunohistochemical staining of Human Li...



Immunohistochemical staining of Human Sm...