



Product Datasheet

IL3RA antibody (orb233662)

Description	Goat polyclonal antibody to IL3RA
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	IL3RA / CD123
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoformsWB: Approx 43kDa band observed in lysates of cell lines Jurkat and Molt4. Recommended concentration: 1-3 µg/ml.
Protein Sequence	RQQYECLHYKTD
Clonality	Polyclonal
MW	43.3; 34.5
Entrez	3563
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat (Lane 1)...



Immunohistochemical staining of Human Co...