


Product Datasheet

HAVCR1 antibody (orb233669)

Description	Goat polyclonal antibody to KIM1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TIM-1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	The immunizing peptide represents part of the extracellular domain. Reported variants represent identical proteinWB: Approx 52kDa band observed in Human Kidney, Testis and Uterus lysates
Protein Sequence	THVTYRKDTRYK
Clonality	Polyclonal
MW	39.2
Entrez	26762
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Kidney (Lane 1)...