

Product Datasheet

Vasopressin peptide (orb233747)

Description	Vasopressin [Arg8] peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Main metabolite of [Arg8]-Vasopressin. It improved the scopolamine-induced deficit of spatial memory in rats about 10'000 fold more than [Arg8]-Vasopressin and therefore has the most potent activity of all the endogenous metabolites of [Arg8]-Vasopressin.
Protein Sequence	Gln-Asn-Cys-Pro-Arg-Gly-NH ₂
Formula	C ₂₅ H ₄₄ N ₁₂ O ₈ S
Purity	> 95%
Source	Synthetic
MW	672.76
CAS Number	96027-30-4
Expiration Date	12 months from date of receipt.

