

Product Datasheet

Dalargin peptide (orb233785)

Description	Dalargin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Dalargin (DALE) has originally been proposed as an antiulcer agent with cytoprotective properties, but it has been tested for other activities with promising results. Dalargin shows antioxidant, antiischemic, and wound-healing properties.
Protein Sequence	Tyr-D-Ala-Gly-Phe-Leu-Arg
Formula	$C_{35}H_{51}N_9O_8$
Purity	> 95%
Source	Synthetic
MW	725.85
CAS Number	81733-79-1
Expiration Date	12 months from date of receipt.

