

Product Datasheet

Enkephalin peptide (orb233799)

Description	Enkephalin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Gly-Gly-Phe-Met-Arg-Arg
Formula	$C_{39}H_{59}N_{13}O_9S$
Purity	> 95%
Source	Synthetic
MW	886.05
CAS Number	76496-10-1
Expiration Date	12 months from date of receipt.

