

Product Datasheet

Oxytocin peptide (orb233802)

Description	Oxytocin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	TGOT exhibits a 640-fold increase in oxytocic/antidiuretic selectivity relative to oxytocin.
Protein Sequence	Cys-Tyr-Ile-Thr-Asn-Cys-Gly-Leu-Gly-NH ₂ (Disulfide bridge Cys1-Cys6)
Formula	C ₃₉ H ₆₁ N ₁₁ O ₁₂ S ₂
Purity	> 95%
Source	Synthetic
MW	940.10
CAS Number	60786-59-6
Expiration Date	12 months from date of receipt.

