

Product Datasheet

Enkephalin peptide (orb233809)

Description	Enkephalin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This Met-enkephalin-containing polypeptide has been isolated from extracts of bovine adrenal medulla.
Protein Sequence	Tyr-Gly-Gly-Phe-Met-Arg-Gly-Leu
Formula	$C_{41}H_{61}N_{11}O_{10}S$
Purity	> 95%
Source	Synthetic
MW	900.08
CAS Number	80501-44-6
Expiration Date	12 months from date of receipt.

