

Product Datasheet

Bradykinin [Leu8, Des-Arg9] peptide [Out of stock] (orb233810)

Description	Bradykinin [Leu8, Des-Arg9] peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	B1 bradykinin receptor antagonist.
Protein Sequence	Arg-Pro-Pro-Gly-Phe-Ser-Pro-Leu
Formula	$C_{41}H_{63}N_{11}O_{10}$
Purity	> 95%
Source	Synthetic
MW	870.03
CAS Number	64695-06-3
Expiration Date	12 months from date of receipt.

