

Product Datasheet

Angiotensinogen peptide (orb233817)

Description	Angiotensinogen peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	As a noncompetitive angiotensin II antagonist, saralasin inhibits the pressor effects of angiotensin II in rats and shows blood pressure lowering activity in humans. It is used as a diagnostic for renin-dependent hypertension.
Protein Sequence	Sar-Arg-Val-Tyr-Val-His-Pro-Ala
Formula	$C_{42}H_{65}N_{13}O_{10}$
Purity	> 95%
Source	Synthetic
MW	912.10
CAS Number	34273-10-4
Expiration Date	12 months from date of receipt.

