

Product Datasheet

Human, Mouse [Des-Asp187, Met186] Melanocyte Protein PMEL 17 peptide [Out of stock] (orb233825)

Description	Human, Mouse [Des-Asp187, Met186] Melanocyte Protein PMEL 17 peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Met-Gln-Val-Pro-Phe-Ser-Val
Formula	C ₄₃ H ₆₉ N ₉ O ₁₁ S
Purity	> 95%
Source	Synthetic
MW	920.15
Expiration Date	12 months from date of receipt.

