

Product Datasheet

Angiotensinogen peptide (orb233826)

Description	Angiotensinogen peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	[Sar1, Ile4, 8]-Angiotensin II activates extracellular signal-regulated kinases 1 and 2 (ERK1/2) via angiotensin II type 1A receptors. Activation is dependent on β-arrestin 2 expression.
Protein Sequence	Sar-Arg-Val-Ile-Ile-His-Pro-Ile
Formula	C ₄₃ H ₇₅ N ₁₃ O ₉
Purity	> 95%
Source	Synthetic
MW	918.15
CAS Number	185461-45-4
Expiration Date	12 months from date of receipt.

