

Product Datasheet

Neurokinin B peptide (orb233902)

Description	Neurokinin B peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This analog is a selective and competitive antagonist of neurokinin B in the guinea-pig ileum.
Protein Sequence	Asp-D-Pro-His-Asp-Phe-D-Trp-Val-D-Trp-Leu-Nle-NH ₂
Formula	C ₆₇ H ₈₇ N ₁₅ O ₁₄
Purity	> 95%
Source	Synthetic
MW	1326.53
CAS Number	109212-72-8
Expiration Date	12 months from date of receipt.

