

Product Datasheet

Mouse, Rat Tyr-CGRP peptide (orb233933)

Description	Mouse, Rat Tyr-CGRP peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	For radioimmunoassay.
Protein Sequence	Tyr-Val-Lys-Asp-Asn-Phe-Val-Pro-Thr-Asn-Val-Gly-Ser-Glu-Ala-Phe-NH ₂
Formula	C ₈₂ H ₁₂₀ N ₂₀ O ₂₅
Purity	> 95%
Source	Synthetic
MW	1785.99
CAS Number	198277-54-2
Expiration Date	12 months from date of receipt.

