

Product Datasheet

CD86 Antibody (orb234297)

Description
CD86 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

E.coli-derived human CD86 recombinant protein (Position: A24-F329). Human CD86 shares 50% amino acid (aa) sequence identity with mouse CD86.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for CD86 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


37682 MW

Uniprot ID
P42081
Dilution Range


Western blot, 0.1-0.5µg/ml Immunohistochemistry(Paraffin-embedded Section), 2-5µg/ml

Expiration Date


12 months from date of receipt.



WB analysis
of CD86 using
anti-CD86
anti...



WB analysis
of CD86 using
anti-CD86
anti...



IHC analysis
of CD86 using
anti-CD86
ant...