

Product Datasheet

NUP98 Antibody (orb234348)

Description

NUP98 Antibody

130KD — 1 2 3 4 5

Species/Host

Rabbit

Reactivity

Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Western blot
analysis
using Anti-
NUP98 P...

Immunogen

E.coli-derived human NUP98 recombinant protein (Position: H549-F880). Human NUP98 shares 95% and 94% amino acid (aa) sequence identity with mouse and rat NUP98, respectively.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for NUP98 is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

197579 MW

Uniprot ID

P52948

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Monkey, Mouse, Rat

Expiration Date

12 months from date of receipt.