

Product Datasheet

PER2 Antibody (orb234358)

Description

PER2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

E.coli-derived human PER2 recombinant protein (Position: N13-K330). Human PER2 shares 84% and 83% amino acid (aa) sequence identity with mouse and rat PER2, respectively.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for PER2 is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

136579 MW

Uniprot ID


O15055

Dilution Range

Western blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



WB analysis
of PER2 using
anti-PER2
anti...