

Product Datasheet

Collagen I peptide (orb234441)

Description	Collagen I peptide
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	COL1A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234810 (Collagen 1 alpha 1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P02452, P11087, P02454
Entrez	1277, 12842, 29393
Expiration Date	12 months from date of receipt.