

Product Datasheet

Collagen V alpha 2 peptide (orb234444)

Description	Collagen V alpha 2 Peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	COL5A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234812 (Collagen 5 alpha 2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05997, Q3U962
Entrez	12832, 1290
Expiration Date	12 months from date of receipt.