

Product Datasheet

DnaJC3 peptide (orb234456)

Description	The peptide is used to anti-DnaJC3 antibody reactivity
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	DNAJC3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234825 (DnaJC3 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q91YW3, Q13217, Q9R0T3
Entrez	63880, 5611, 100037258
Expiration Date	12 months from date of receipt.