

Product Datasheet

DNMT3A peptide (orb234457)

Description	The peptide is used to anti-DNMT3A antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DNMT3A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234826 (DNMT3A antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q1LZ53 , Q9Y6K1 , O88508
Entrez	13435 , 1788 , 444984
Expiration Date	12 months from date of receipt.