

Product Datasheet

EDEM1 peptide (orb234459)

Description	The peptide is used to anti-EDEM1 antibody reactivity
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	EDEM1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234828 (EDEM1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92611 , Q925U4
Entrez	192193 , 9695
Expiration Date	12 months from date of receipt.