

Product Datasheet

HLA-A Blocking Peptide (orb234503)

Description	HLA-A Blocking Peptide
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	HLA-A; HLA-B; HLA-C
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234872 (HLA-A/B/C antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P30456, P30512, P10314, P01892, P01891, P04439, P18462, P30443, P30453, Q09160, P10316, P16189, P30455, P05534, P30450, P30457, P13746
Entrez	3106, 3107, 3105
Expiration Date	12 months from date of receipt.