

Product Datasheet

TSHZ3 peptide (orb234606)

Description	The peptide is used to anti-TSHZ3 antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TSHZ3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234975 (TSHZ3 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63HK5, Q8CGV9, D3ZKB9
Entrez	243931, 308523, 57616
Expiration Date	12 months from date of receipt.