

Product Datasheet

ADAM28 peptide (orb234626)

Description	The peptide is used to anti-ADAM28 antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	ADAM28
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb234995 (ADAM28 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JLN6, Q9UKQ2
Entrez	13522, 10863
Expiration Date	12 months from date of receipt.