

Product Datasheet

ALDH1A1 peptide (orb234635)

Description	ALDH1A1 Peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ALDH1A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb235004 (ALDH1A1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P24549, P00352, P51647
Entrez	11668, 24188, 216
Expiration Date	12 months from date of receipt.